

Spectrolaser Application Note – Steel Industry


The Spectrolaser's use of Laser Induced Breakdown Spectroscopy (LIBS) is well known in its ability to achieve accurate and fast elemental analysis of materials.

In the steel industry the properties of LIBS that are particularly beneficial include:

- The ability to analyze conductive and non-conductive materials in the one instrument.
- The ability to analyze small samples and inclusions in samples using the focused laser beam,
- The ability to “burn off” crusts and surface impurities using the ablating power of the laser prior to analysis, and
- The ability to obtain elemental mapping across a surface using XY sample scanning.

A broad range of elements can be accurately measured in steel and slags using the Spectrolaser including both major components (Fe, C, Cr, Ni, Ti...) and minor components (K, Ca, Mg, Al, Si ...).


For slag samples the material is typically crushed and pressed into pellets prior to analysis.



Being of compact size and requiring no services apart from power, the Spectrolaser can easily be located close to the process thereby improving sample analysis times and lowering production costs.

Scanning and Imaging Applications

With the Spectrolaser Target's scanning and imaging ability there are many further applications. In the example below variations in steel plate are routinely analyzed with a Spectrolaser Target on a factory floor.



Analysis Benefits

- a low level of sample preparation is required improving lab efficiency and lowering production costs,
- simultaneous analysis of all detectable elements,
- analyze conductive and non-conductive materials in the one instrument.
- analyze small samples and inclusions in samples, and
- analyze coating or alternatively use “burn off” prior to analysis.

Spectrolaser Features

- Intuitive design allows easy operation in any environment,
- Complete software control over laser power, detector exposure and analysis times,
- Lab manager restricted access to instrument settings,
- Average multiple laser pulses,
- Overlay spectra for easy comparison,
- Peak searching using a database of >10,000 emission lines,
- Multiple report saving and export options
- Graphical interface for elemental mapping (Spectrolaser Target XY)
- Fully supported worldwide